

# SiMPLi™ Series

## Ultra Simple USB Ultrasound

Connects to any Windows Computer, Laptop or Tablet  
Super Lightweight, Handheld and Portable  
Plug & Play – Instant Real Time Imaging  
100% Manufactured in the USA

### Images

- Bladder
- Kidney
- Spleen
- Heart
- Liver



### Procedures

- Fine-needle aspiration
- Pregnancy Diagnosis
- Abdominocentesis
- Cystocentesis
- Liver Biopsy


**INTERSON**  
medical instruments

7150 Koll Center Parkway  
Pleasanton, CA 94566, USA

Tel: (925) 462-4948  
Fax: (925) 462-4833  
sales@interson.com  
www.interson.com

# SP-L01 Small Parts Pet Scan

## Specifications



Model Number	99-7912
Depth Range	0.1 cm - 10 cm
Focus	Dynamic over full depth range
Pulse Frequency	5 - 10 MHz
Patient Contact	38 mm
Cable	USB 2.0, 2 meter (5 meter opt.)
Size	142 mm x 58 mm
Weight	80 grams
System Requirements	OS: Windows 7, 8, 10 Processor / RAM: 2 GHz / 4 GB

Interson designs, manufactures, sells, and supports a wide variety of USB ultrasound probes in our ISO 13485 facility. All probes are manufactured in the United States and are subject to regulation by the U.S. Food and Drug Administration and other international, state, and local regulations. Interson's SeeMore Windows application and MS Visual Studio SDK provide a real time display of ultrasound images on desktop, notebook, and tablet computing devices supporting end-users and research. Interson's USB ultrasound is the ideal solution for handheld point-of-care use in human and veterinary applications. Our standards-based WebRTC interface is ideal for telemedicine and remote diagnostic applications. All probes come with a standard two year warranty, extended warranties are available.



7150 Koll Center Parkway  
Pleasanton, CA 94566, USA

Tel: (925) 462-4948  
Fax: (925) 462-4833  
sales@interson.com  
www.interson.com