


# ViewBladder 10

## *Point-of-Care Bladder Volume*



- Long Term Care
- Rehab
- Acute Care
- Hospital Post-Op

## Tablet Bladder Scanner

- ***Ready to go right out of the box !***
- **Post-Void Residual (PVR)** – reduce catheterizations, UTIs
- **Rehab** – 80% of patients can be rehabbed in 45 days or less
- **Hospital Post-Op** – accurately assess retention
- **Real-Time Image** – visually determine full versus empty
- **Accurate Measurements** – 0-1500 ml., +/-10% +/- 10 ml

# Tablet Bladder Scanner

---

Model Number	99-9008
Tablet	Hewlett Packard Pavilion X2 or equivalent
Operating System	Windows 8.1 Home
Tablet Size	17.5 cm x 26.3 cm x 1.6 cm
Tablet Weight	900 grams
Measurement Range	0-1500 ml, (ccm)
Measurement Accuracy	+/- 10% +/- 10 ml
Depth Range	4 - 24 cm
Pulse Frequency	2.5 – 5.0 MHz
Scan Angle	90 degrees
Patient Contact	35 mm
Frame Rate	12 fps
Focal Point	7.5 cm
Probe Size	17.3 cm x 4.1 cm
Probe Weight	250 grams
Cable	USB 2.0, 2 meter (5 meter opt.)

Interson designs, manufactures, sells, and supports a wide variety of USB ultrasound probes in our ISO 13485 facility. All probes are manufactured in the United States and are subject to regulation by the U.S. Food and Drug Administration and other international, state, and local regulations. Interson's SeeMore Windows application and MS Visual Studio SDK provide a real time display of ultrasound images on desktop, notebook, and tablet computing devices supporting human and veterinary applications. Interson's USB ultrasound is the ideal solution for telemedicine, and research. All probes come with a standard one year warranty, extended warranties are available.