


USB Ultrasound Probe

Transvaginal 99-5909



Left Ovary - Sagittal

Right Ovary - Transverse

Applications

- Ectopic Pregnancy
- Valoration – Cervix, Ovaries
- Douglas Pouch – free fluid
- Cyst Presence / Volume

Specifications

Depth Range: 0.5 cm - 10 cm	Focal Zone: 0.5 cm - 4 cm
Pulse Frequency: 7.5 MHz, 12 MHz, 15 MHz, 24 MHz	Focal Point: 2 cm
Scan Angle: 90 degrees	Patient Contact: 21 mm
Frame Rate: 15 fps	Cable: USB 2.0, 2 m (5 m opt.)
Size: 15.5 cm by 3.5 cm	Computer OS: Windows 7, 8
Weight: 100 g	Processor, RAM: 4 GHz, 4 MB

Interson designs, manufactures, sells, and supports a variety of USB ultrasound probes in our ISO 13485 facility. All probes are manufactured in the United States and are subject to regulation by the U. S. Food and Drug Administration and other international, state, and local regulations. Interson's SeeMore Windows application and SDK provide procedure assisted imaging and real time display of ultrasound images on desktop, notebook

and tablet computing devices for human, veterinary, telemedicine, and research applications. Probes come with a standard one year warranty, extended warranties are available.