

SeeMore™ USB Probes



We **revolutionize** imaging devices.

CALCULATION AND MEASUREMENT PACKAGES

OB/GYN Calculation Package (optional add-on)

- **CR** (Crown Rump)
- **GS** (Gestational Sac)
- **FL** (Femur Length)
- **HC** (Head Circumference)
- **AC** (Abdominal Circumference)
- **BPD** (Bi-Parietal Diameter),
- **AFI** (Amniotic Fluid Index)
- Bladder Volume

Urology Calculation Package (optional add-on package)

- Bladder Volume
- Prostate Volume

Standard Measurements and Annotation Package (included with all software)

- Distance
- Area
- Annotation
- Pointers