

# SeeMore™ USB Probes



We **revolutionize** imaging devices.

## GENERAL PURPOSE ABDOMINAL PROBES



### AB 3.5 MHz

- Center Frequency: 3.5 MHz
- Pulse Frequency : 3.5 MHz and 5.0 MHz
- Scan Angle: 90°
- Focal Point: 7.5 cm
- Scan Depth (Max): 20 cm
- Clinical Applications : Abdomen (kidney, liver, gall bladder, abdominal aorta), Uterus, Bladder, FAST exam

### AB 5.0 MHz

- Center Frequency: 5.0 MHz
- Pulse Frequency : 5.0 MHz and 3.5 MHz
- Scan Angle: 90°
- Focal Point: 6 cm
- Scan Depth (Max): 20 cm
- Clinical Applications : Pediatric and smaller patients