



SiMPLi Series

SP-L01 Point-of-Care USB Linear Array Transducer

Ready for user-friendly ultrasound at the patient bedside? Simply connect the USB probe cable directly to your Windows laptop or tablet and you instantly have a lightweight, handheld and portable ultrasound imaging system for all of your point-of-care imaging needs.



Specifications

Depth Range	0.1 cm - 10 cm
Focus	Dynamic over full depth range
Bandwidth	5 - 10 MHz
Features	2D, linear array, color flow imaging, lightweight transducer with flexible cable
Indications	Peripheral Vascular, Small Parts, Pediatrics, MSK, Neuromuscular, Pain Management, Acute Care
Weight	80 grams
Cable	USB 2.0, 2 meter (5 meter optional)
System Requirements	OS: Windows 7, 8, 10 Processor/RAM: 2 GHz/4 GB

Lightweight. Reliable. Flexible.

For more information regarding the SP-L01 or our other transducers, contact us at **interson.com/contact-us** or call us at **(925) 462-4948**.



Brachial Access



Thyroid/Carotid/Jugular



Carotid/Intima Media

