



ViewBladder 10

Point-of-Care Bladder Volume



Ready for user-friendly ultrasound at the patient bedside? Simply connect the USB probe cable directly to your Windows laptop or tablet and you instantly have a lightweight, handheld and portable ultrasound imaging system for your point-of-care imaging needs.

Interson is the leader in USB ultrasound technology. From here came the development of ViewBladder 10, specifically designed for assessing bladder volume. A simple solution. Bladder scanning at your fingertips.



Specifications

Probe Model Name	GP 3.5 USB Ultrasound Probe
Part Number P/N	99-5929
Measurement Range	10-1000 ml, (ccm)
Depth Range	4-24 cm
Pulse Frequency	2.5-5.0 MHz
Scan Angle	90 degrees
Patient Contact	35 mm
Frame Rate	12 fps
Probe Size	17.3 cm x 4.1 cm
Probe Weight	250 grams
Cable	USB 2.0, 2 meter (5 meter opt.)

Tablet Specifications

Tablet	Intel/AMD X86_64 - based tablet
Operating System	Windows 10 or newer (64 bit only)
Tablet Size	10" Screen diagonal or similar
Screen Features	Touch screen preferred
Processor	Quad Core CPU or better

Affordable. Reliable. Flexible.

For more information regarding the VB10 or our other transducers, contact us at interson.com/contact-us or call us at (925)462-4948.

- **Post-Void Residual (PVR)** – Ability to view and measure bladder volume below 50 ml
- **Rehabilitation** – Significant improvement in the reduction of rehab time
- **Hospital Post-Op** – Accurately assess retention
- **Real-Time Image** – Visually determine bladder size and shape

Low cost of ownership. No need for yearly calibration. Affordable service fee. 1 year limited warranty.

Ready to go right out of the box.

ViewBladder 10 Kit (P/N 99-9008)

Probe and Tablet included with pre-loaded VB10 app.

